

UNIVERSITY OF CALIFORNIA,
IRVINE

Spin torque driven magnetization dynamics in nanoscale magnetic tunnel junctions

DISSERTATION

submitted in partial satisfaction of the requirements
for the degree of

DOCTOR OF PHILOSOPHY

in Physics

by

Chengcen Sha

Dissertation Committee:
Professor Ilya Krivorotov, Chair
Professor Zuzanna Siwy
Professor Wilson Ho

2018

© 2018 Chengcen Sha

DEDICATION

To My parents, Zhenglian and Wenyu.

TABLE OF CONTENTS

	Page
LIST OF FIGURES	iv
LIST OF TABLES	v
ACKNOWLEDGMENTS	vi
ABSTRACT OF THE DISSERTATION	vii
Bibliography	1
A Appendix Title	2
A.1 Detailed System Design of Perpendicular station	2

LIST OF FIGURES

Page

LIST OF TABLES

Page

ACKNOWLEDGMENTS

I would like to express my deepest appreciation to Professor Ilya Krivorotov. It is quite a privilege to work with him. I really learned a lot from his great knowledge of physics and sharp insights into problem solving. I still remember almost six years ago I had my first Skype interview with Ilya and he asked me about exchange bias effect. Without his guidance and help this dissertation would not be possible.

I would also like to thank all my lab mates : Zheng Duan, Igor Barsukov, Eric Montoya, Brian Youngblood, Yu-Jin Chen, Andrew Smith, Han Kyu Lee, Jenru Chen, Jieyi Zhang, Alejandro Jara, Chris Safranski and Josh Dill. I have been received a great amount of help since the first day I joined the lab. When I first joined the lab, I did not know even most simple thing, such as using a torque wrench and basis soldering. These guys helped me overcome a lot of difficulties and provided help whenever I asked. I never wished I could be surrounded by such a amazing group of people and I really feel honored to call these guys my lab mate.

During the last five years in the United States, I am happy to keep in touch with my good friends over the world. I would like to thank Jiao Li, Rui Da, Ruohui Yang, Xun Liu, Hongyu Zhu and Qingyu Zhu. Those people are all pursuing or have already obtained their Ph.D. degree and we have had lots of useful/useless discussions about our common experiences. We shared numerous joy and distress over the past five years. I might not be able to see these friends quite often and in fact, I have not met some of them for almost five years. However, their support is greatly valuable to me and I really hope nothing but the best for all of them from my heart.

Of course I am extremely luck to have a loving family. My parents have been the greatest support to me without any conditions. They have provided all they can give to me and I cannot pay back them enough. My wife has been the angle of my life and I love her very much. Of course I cannot forget my cat Luca, who gives me lots of joy when I was lonely.

ABSTRACT OF THE DISSERTATION

Spin torque driven magnetization dynamics in nanoscale magnetic tunnel junctions

By

Chengcen Sha

Doctor of Philosophy in Physics

University of California, Irvine, 2018

Professor Ilya Krivorotov, Chair

Change the abstarcet!!!!!!

Bibliography

Appendix A

Appendix Title

A.1 Detailed System Design of Perpendicular station

We have developed one of the most sensitive Spin-torque ferromagnetic resonance perpendicular magnetic stations. Here we describe the detailed system designs.

First of all, here is all the equipments needed to build the out-of-plane magnetic probe station:

1. GMW Dipole Electromagnet Model 3470
2. Kepco bipolar operational power supply model Model 50-8M
3. Cascade RF probe : SG-120um
4. Cascade RPP210-AI probe positioner (both the probe and the positioner are non-magnetic)
5. Sentech Output 720p Cased Camera

Perpendicular ST-FMR station




Figure A.1: Perpendicular ST-FMR station Setup

6. Navitar 12X Zoom Lens System

7. AmScope LED-80M 80-LED Microscope Ring Light

Fig.A.1 sketches the design of the out-of-plane station. The magnet is fixed vertically on metal frame and the stage height is adjustable. At first, we can land the probe to make contact of the sample with the magnet moving away(shown in Fig.A.2(a)). After making contact of the sample, first remove the camera (the setup could be improved by making a stationary camera). Then we can slide in the magnet so that the sample is located in the center of the magnet. It is important not to touch the probe and microwave cable when sliding the magnet. After moving the magnet we are ready to make ST-FMR measurement as showing in Fig.A.2(b).

Here is the list of equipments needed for making ST-FMR measurements up to 40 GHz.




Figure A.2: Operation of Out-of-plane probe station




Figure A.3: (a) Side view of the magnet. (b) Bottom view of the magnet

1. Signal recovery 7225 DSP lock-in Amplifier
2. Hittite HMC-T2240 synthesized signal generator, 10 MHz to 40 GHz
3. Clear Microwave Broadband Bias Tee BT50K40 50Khz-40GHz
4. Keithley 2400 Source Meter (DC Source)
5. Microwave cable : Teledyne Accutest R95-0004-072 (72 inch 1GHz- 40GHz)
6. Pomona BNC cables

When connecting the microwave cables, there are two small things worth notice. Firstly, all the connectors should be cleaned regularly. Secondly, the microwave cables should be supported to release all the possible tensions.

The last pieces of equipments needed is the field modulation setup. Here is the list:

1. Behringer EUROPOWER Professional 4,000-Watt Stereo Power Amplifier
2. TE CONNECTIVITY / CGS CJT10004R7JJ Through Hole Wire wound Resistor, 4.7 Ohm
3. High quality cable connecting from lock-in to the input of the sound amplifier : Monster Performer 500 - 10' Speaker Cable
4. BK Precision 2831E Ammeter to control the current through the copper wire

When making the field-modulation coil, it is better to use any low resistance copper wire for magnetic field modulation coil. For out-of-plane field modulations, we use an external coil above the sample as shown in Fig.A.4(a). In-plane field modulation can be achieved by simply placing straight wire above the sample as shown in Fig.A.4(b). In our current design,




Figure A.4: (a) Coil for out-of-plane modulation (b) Wire for in-plane modulation

we embed the coil into the plastic base right under the sample. It is important to ensure that the wire is not in contact with the sample or any other parts of the microwave setup.

In the experiment, two parameters need to be determined for the field modulation: the amplitude of modulation field and the frequency of the modulation ac field. Typically Increasing the modulation field increases the signal amplitude, but not infinitely as shown in Fig.A.5(a). If you are over-modulating, the signal becomes distorted and the linewidth broadens as shown in Fig.A.5(b). In our current setup, the input ac current is about 3.6 A to achieve the modulation field around a few oersted field, which is enough to have decent signal-to-noise ratio without distorting the spectrum.




Figure A.5: (a) Peak-to-peak amplitude versus rms modulation field over real peak-to-peak linewidth. (b) peak-to-peak linewidth over real peak-to-peak linewidth versus rms modulation field over real peak-to-peak linewidth.

Mode 0	10.98	11.92	11.57	11.98	10.94	11.71	10.94	11.84	11.38	11.8	12.09
Mode 1	12.85	14.6	13.98	14.64	13.21	14.28	13.05	13.6	14.29	14.66	14.01
Mode 2	13.75	15	15.4	15.33	13.2	15.48	13.66	14.57	14.86	15.24	14.18
Gap 0-1	1.87	2.68	2.41	2.66	2.27	2.57	2.11	1.76	2.91	2.86	1.92
Gap 0-2	2.77	3.08	3.83	3.35	2.26	3.77	2.72	2.73	3.48	3.44	2.09
Gap (1+2)/2 - 0	2.32	2.88	3.12	3.005	2.265	3.17	2.415	2.245	3.195	3.15	2.005

Table A.1: C06 70nm summary

Mode 0	9.83	10.5	11.05	10.54	11.89	11.39	11.2	10.39	10.61	10.39	11.07
Mode 1	11.76	12.68	13.24	12.61	13.78	13.6	13.19	12.85	13.47	12.6	13.02
Mode 2	12.71	13.04	13.86	13.04	14.88	14.29	13.5	14.06	13.79	12.99	13.99
Gap 0-1	1.93	2.18	2.19	2.07	1.89	2.21	1.99	2.46	2.86	2.21	1.95
Gap 0-2	2.88	2.54	2.81	2.5	2.99	2.9	2.3	3.67	3.18	2.6	2.92
Gap (1+2)/2 - 0	2.405	2.36	2.5	2.285	2.44	2.555	2.145	3.065	3.02	2.405	2.435

Table A.2: C07 80nm Summary

Mode 0	9.94	10.19	9.28	10.54	9.14	9.14	9.6	9.77	9.13	10.51
Mode 1	11.65	12.41	10.93	12.29	11.2	11.2	12.07	11.8	10.69	11.96
Mode 2	12.15	12.66	12.07	12.29	12.05	12.05	12.36	12.9	11.26	12.65
Gap 0-1	1.71	2.22	1.65	1.75	2.06	2.06	2.47	2.03	1.56	1.45
Gap 0-2	2.21	2.47	2.79	1.75	2.91	2.91	2.76	3.13	2.13	2.14
Gap 0-(1,2)	1.96	2.345	2.22	1.75	2.485	2.485	2.615	2.58	1.845	1.795

Table A.3: C08 90nm summary

Mode 0	7.88	9.62	9.09	8.97	9.07	8.84	9.24	8.35	8.75	8.77	7.96
Mode 1	9.55	11.4	10.43	10.36	10.4	10.25	10.47	10.23	10.35	10.57	9.89
Mode 2	9.96	11.59	10.52	11.2	10.7	10.52	11.8	10.59	11.31	10.61	10.49
Gap 0-1	1.67	1.78	1.34	1.39	1.33	1.41	1.23	1.88	1.6	1.8	1.93
Gap 0-2	2.08	1.97	1.43	2.23	1.63	1.68	2.56	2.24	2.56	1.84	2.53
Gap 0-(1,2)	1.875	1.875	1.385	1.81	1.48	1.545	1.895	2.06	2.08	1.82	2.23

Table A.4: C09 120 nm

Mode 0	8.02	7.61	7.81	8.13	7.86	7.94	7.85	7.46	7.18	8.33	8.73	7.98	7.65
Mode 1	9.71	8.8	9.36	9.13	9.27	9.37	9.16	8.79	8.13	9.3	9.73	9.14	8.77
Mode 2	9.9	8.96	9.83	9.81	9.52	9.37	9.44	9.42	8.71	9.91	10.22	9.6	9.21
Gap 0-1	1.69	1.19	1.55	1	1.41	1.43	1.31	1.33	0.95	0.97	1	1.16	1.12
Gap 0-2	1.88	1.35	2.02	1.68	1.66	1.43	1.59	1.96	1.53	1.58	1.49	1.62	1.56
Gap 0-(1,2)	1.785	1.27	1.785	1.34	1.535	1.43	1.45	1.645	1.24	1.275	1.245	1.39	1.34

Table A.5: C10 150 nm summary